



MATERIALS:  
 TOP PLATE - 6061 T6 ALUMINUM  
 BASE PLATE - STAINLESS STEEL  
 MAGNETS - NDFEB

NOTE: FINISH TO BE CLEAR HC PER SAMPLE

PROPRIETARY AND CONFIDENTIAL	REV.	DATE	NAME	PERMAGEN LABWARE	
THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF PERMAGEN LABWARE, LTD. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF PERMAGEN LABWARE IS PROHIBITED.	A	9.17.16	KAB	TITLE: 24-WELL MAG PLATE	
				DWG.NO.: 24001	
				MATERIAL:	REV.: A
				SCALE: 1 : 1	FINISH: SEE NOTE