



NOTES:
MATERIALS - ALUMINUM, STAINLESS STEEL, NdFeB
MAGNETS
FINISH - CLEAR HC ANODIZE

PROPRIETARY AND CONFIDENTIAL	REV.	DATE	NAME	PERMAGEN LABWARE	
THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF PERMAGEN LABWARE, LTD. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF PERMAGEN LABWARE IS PROHIBITED.				TITLE: P384 POST MAGNET PLATE 384	
				DWG.NO.: A071538	
				MATERIAL: SEE NOTE	REV.:
				SCALE:	FINISH: SEE NOTE