

# 96S Super Magnet Plate 96-Well




Permagen's 96S Super Magnet plate was designed as a direct replacement magnet for the more expensive 96S plates on the market. Used for manual or automated applications from most common microplates. Contains strong NdFeB ring magnets for separations up to 2 mL, or faster small volumes. Includes integrated spring cushion base for improved recovery by compensating for robot, and labware inconsistencies.

## Features

- Guaranteed to meet or exceed performance of either Part No.: A32782 or A001322\*
- Almost 50% in cost savings\*\*
- Spring Cushion base compensates for labware and robot inconsistencies for precise pipetting and maximum recovery
- ANSI/ SBS Footprint to fit standard liquid handler nests

\*Uses same magnets, features all of the same dimensions, and same materials  
\*\*When comparing our MSRP of \$585.00 to Beckman® \$1131.00



## Compatible...

Full-, Semi-, and Non-Skirted PCR Plates  
Deep-Well Plates  
Standard Plates  
Round-, and V-Bottom Plates



Products can be ordered worldwide at our online shop, by phone, or email

Cat No.	Product	Product Name
A092722	96-Well Ring Magnet Plate	96S Super Magnet



Permagen Labware, LTD  
110 Haverhill RD. STE 320  
Amesbury, MA 01913  
P: +1 978-914-4943

# X96 Dual Ring Magnet Plate




Permagen Labware, LTD  
110 Haverhill RD. STE 320  
Amesbury, MA 01913  
P: +1 978-914-4943