




R 1.52

85.47 (SBS/ANSI)

Ø 6.15 ±.05mm

127.25 (SBS/ANSI)

110.32

127.76 (SBS/ANSI)

74.32

14.00

10.01

PROPRIETARY AND CONFIDENTIAL	REV.	DATE	NAME	PERMAGEN LABWARE	
THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF PERMAGEN LABWARE, LTD. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF PERMAGEN LABWARE IS PROHIBITED.				TITLE: 96-WELL CHILLER	
				DWG.NO.: 01961	
				MATERIAL: 6061	REV.: B
				SCALE:	FINISH: SEE NOTE

DIMENSIONS ARE IN MM

NOTES: PLATING TO BE GRAY HC ANODIZED PER SAMPLE